

## Making LB Agar Plates

Batch makes about 40 plates.

### **Making the LB Agar**

1. Add 250 mL of dH<sub>2</sub>O to a graduated cylinder.
2. Weigh out 20g of premix LB Agar powder (VWR DF0445-17) or:
3. 5.0 g tryptone
4. 2.5 g yeast extract
5. 5.0 g NaCl
6. 7.5 g agar
7. Mix powder well to bring into solution
8. Add dH<sub>2</sub>O to total volume of 500 mL and transfer to 1 L flask
9. Put on stirring hot plate and heat to boil for 1 min while stirring.
10. Transfer to 1 L pyrex jar and label with autoclave tape.
11. Autoclave at liquid setting for 20 minutes in a basin making sure to loosen top
12. Let agar cool to ~55C (you should be able to pick up the jar without a glove)

### **Pouring the Plates**

1. Make sure bench top has wiped down with bleach/EtOH.
2. Remove sterile Petri dishes (VWR 25384-208) from plastic bag (save the bag for storage).
3. Pour a thin layer (5mm) of LB Agar (~10mL) into each plate being careful to not lift the cover off excessively (you should be able to just open up enough to pour).
4. Swirl plate in a circular motion to distribute agar on bottom completely.
5. Let each plate cool until its solid (~20 minutes) then flip so as to avoid condensation on the agar.
6. Store plates in plastic bags in fridge with: name, date and contents (note any additive).

### **Special Additives** (to be added to LB Agar right before pouring plates)

Ampicillin (VWR 80055-786)      50 mg dissolved in a small amount of dH<sub>2</sub>O  
(concentration 100 ug/mL)

X-gal (VWR IB02260)      50 mg dissolved in a small amount of DMSO

### **Stock solutions**

Ampicillin    20mg/mL      200mg in 10mL dH<sub>2</sub>O (store at 4 in 1mL aliquots) use  
50uL on each plate

IPTG (VWR EM-5800)      100mM      238 mg IPTG in 10mL dH<sub>2</sub>O (store at -20  
in 1mL aliquots) use 40uL on each plate

X-gal    40 mg/mL      400 mg X-gal in 10mL DMSO (store at -20 in 1mL aliquots foil  
wrapped tubes) use 40uL on each plate